



# 6702 Hercules

Surgical Table

# Giving operating rooms the ability to do more.



Combining impressive 1,200 pound lift capacity with exclusive table top rotation, the 6702 Hercules exceeds all standards in general purpose surgical tables. Skytron offers 210° table top rotation – making imaging and patient positioning easy and safe. From pediatric to bariatric, the 6702 Hercules gives any operating room the versatility to handle more procedures and specialties.

## Imaging Efficiency

At the turn of a handle, the 6702 Hercules allows superior imaging access at both the head and foot end of the table. The table can quickly rotate toward the C-arm and the low-profile base stays clear of the imaging equipment. Skytron takes the work out of imaging setup.

## Surgical Positioning Flexibility

A remarkable 1,200 pound lift and 1,000 pound articulation capacity makes the 6702 Hercules ready for patients of diverse sizes. A low 23<sup>5</sup>/<sub>8</sub>" table height provides great access to larger patients. The removable back and leg sections allow seamless integration with accessories such as the power beach chair positioner – providing automatic, one-touch setup for shoulder surgery.

## Improved Patient Safety

With a table top that rotates independent of the base, the 6702 Hercules allows the brakes to remain firmly in place while moving the patient and maximizing surgical access. This unique feature sets the 6702 Hercules apart and makes any operating room safe.



Increased  
Surgical  
Specialties



Generous  
Weight  
Capacity



Efficient  
Top  
Rotation



## 6702 Hercules Series

### 210° Top Rotation, General Purpose Table

Skytron offers true top rotation, allowing the surgeon or C-arm to remain at the foot end of the table where the base does not interfere. And with the removable back and leg section, the power beach chair can be attached for the most comfortable shoulder surgery setup available.



### Differentiating Features

- 210° top rotation
- 1,200 lb lift and 1,000 lb articulation weight capacity
- Removable back and leg sections for quick accessory attachment and positioning
- Power beach chair positioner for one-touch set-up
- Durable stainless steel base is easy to clean
- 23 5/8"-41 1/2" adjustable table height
- 30° lateral tilt
- 30° Trendelenburg and reverse Trendelenburg



## One Table, More Specialties

### Orthopedics – Shoulders

- 180° top rotation - position closer to surgical site
- Removable back section - allows attachment for table powered beach chair

### ENT

- No table base under patient - easier access
- Larger area under patient - perform more surgeries while seated

### Imaging

- Top rotation - allows more space for C-arm
- Carbon fiber extension accessory - completely radiolucent from side to side

### Bariatric

- 1,200lb lift - exceeds most bariatric needs
- 1,000lb articulating capacity - easily obtain patient position

### Urology

- 23" of imaging length - simplifies imaging of the kidneys, ureter and bladder
- Light weight transfer board accessories - easy for staff to maneuver



23<sup>5/8</sup>" Minimum  
Table Height



41<sup>1/2</sup>" Maximum  
Table Height



40° Flex



90° Reflex



20° Leg Up



100° Leg Down



One-touch positioning functions include return-to-level and beach chair



Battery life/diagnostic indicator



Easy to use clip style power cord



Optional side extensions widen the table top to 30"



Optional 40" carbon fiber leg section for upper and lower body imaging



Self-leveling brakes for stability on uneven floors

## One-Touch Pendant Control

- Simple one-touch operation
- Backlit control for high visibility in low lighting
- Anatomically correct return-to-level
- Anatomically correct beach chair position
- Recessed buttons prevent accidental activation
- Durable and water resistant
- Can be used in conjunction with optional three function foot control



90° Back Up



40° Back Down



30° Trendelenburg



30° Reverse Trendelenburg



30° Tilt - Left



30° Tilt - Right

# Positioning Guide



Upper Body Imaging



Lower Body Imaging



Lap Nissen



Urology



Ophthalmic/ENT



Nephrectomy

# Bariatric Options



Abdominal



Gall Bladder



Neuro (Neck)



Neuro (Lumbar)



Shoulder Arthroscopy



Neuro (Head)



Hip Revision



Cardiovascular



Cysto/Gyn

For a full offering of Table Accessories please visit [skytron.us/products/pages/table\\_accessories](http://skytron.us/products/pages/table_accessories)



# Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings  
Washers, Ultrasonics, and Cart Washers  
Prep & Pack Tables & Sterilizers



## Clinical

Surgical Tables  
Booms & Lights  
Full Offering of Stainless Products



## Architectural & Environmental

Modular Walls  
Custom Window and Door Offering  
UVC Disinfection Products



## Clinical Business Intelligence

Video Integration  
Asset Tracking & Supply Chain Solutions  
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions. **To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV3 09-22  
PN L2-010-77