


DOVERA SURGICAL LIGHTING

PRODUCT BROCHURE



Suited for any surgical application, DoVera provides superior illumination while eliminating glare.





Double Reflective/ No glare

DoVera's unique Double Reflection technology eliminates glare at the lighthead and provides uniform light across the surgical field. High intensity LED light beams project onto a segmented mirror, which reflects the light onto a smooth second mirror and subsequently to the surgical field. Double Reflection eliminates glare produced by traditional direct lighting technologies



Personalised Lighting

Spot size and Colour Temperature- Capacitive Touch control at the yoke allow for easy, personalised set up of the light. Three spot sizes provide the right amount of light for every procedure while the seven colour temperature options create a range of options from soft to bright white light. The surgical team can adjust the intensity to their preference throughout the procedure via the sterile centre handle. Optional wall control available for adjustments outside of the field.

DoVera Surgical Lighting 👌 Hipac



DoVera Surgical Lighting Features



Intuitive controls on the yoke provide quick personalisation of spot size, colour temperature, and intensity.

Safety Light Safety lighting for darkened surgical suites



Optional Wall Control Non-sterile option for lighting adjustments



Single Yoke Easy positioning for low ceiling heights



Double Yoke Maximum positioning versatility

