

Gel Pressure Care Products

GEL PRESSURE CARE PRODUCTS

The range of Hipac gel products are used in thousands of healthcare facilities around the world. They are manufactured from a high quality polymer gel that is tissue equivalent and emulates human skin, preventing patients from decubitus ulcers and pressure injuries.



Benefits of Hipac gel products include:

- Latex free for patient and staff safety.
- Leak proof – even if punctured. Hipac's advanced polymer gel does not run and as a result has a superior lifespan.
- Impervious to body fluids and microbial growth.
- Easy to clean. Can be wiped down and disinfected with any hospital approved disinfectant.
- Won't bottom out.
- Reusable.

GEL 3 PIECE TABLE PAD SET

- Set contains one head section, one mid section and one foot section, all without cutout.

CODE:



111729 3 Piece Table Pad Set, Gel, no Cutout

GEL 3 PIECE TABLE PAD SET WITH CUTOUT

- Set contains one head section, one mid section with cutout and one foot section with cutout.

CODE:

137893 3 Piece Table Pad Set, Gel, with Cutout



GEL 3/4 TABLE PAD

CODE:

111647 3/4 Table Pad, Gel, no Cutout,
1170[L] x 505[W] x 13mm[H]



GEL 3/4 TABLE PAD WITH CUTOUT

CODE:

111646 3/4 Table Pad, Gel, with Cutout,
1170[L] x 505[W] x 13mm[H]



GEL FULL TABLE PAD

CODE:

111645 Full Table Pad, Gel,
1830[L] x 505[W] x
13mm[H]



GEL TABLE PAD FOOT SECTION

CODE:

111648 Table Pad Foot Section,
Gel, 505[L] x 505[W] x
13mm[H]



GEL TABLE PAD FOOT SECTION WITH CUTOUT

CODE:

111644 Table Pad Foot Section,
Gel, cutout, 508[L] x
505[W] x 13mm[H]



Gel Pressure Care Products

GEL TABLE PAD TORSO SECTION

CODE:


111731 Table Pad Torso Section, Gel, 1000[L] x 505[W] x 13mm[H]



GEL TABLE PAD HEAD SECTION

CODE:

111730 Table Pad Head Section, Gel, 230[L] x 505[W] x 13mm[H]



GEL CONCAVE ARMBOARD PAD

- Contoured for multi direction support.

CODE:

111495 Concave Armboard Pad, Gel, 530[L] x 150[W] x 50mm[H], Large
111496 Concave Armboard Pad, Gel, 410[L] x 150[W] x 50mm[H], Small



GEL DURABOARD PAD

- Use on standard arm board extensions.
- Protects hand, forearm and elbow.

CODE:

111488 Duraboard Pad, Gel, 660[L] x 152[W] x 10mm[H]
111489 Duraboard Pad, Gel, 380[L] x 150[W] x 10mm[H]



GEL ARMBOARD PAD

- Use on standard armboard extensions.
- Protects hand, forearm and elbow.

CODE:

111486 Armboard Pad, Gel, 380[L] x 100[W] x 6mm[H], Small
111487 Armboard Pad, Gel, 635[L] x 125[W] x 6mm[H], Large
139820 Armboard Pad, Gel, 600[L] x 200[W] x 13mm[H], to suit Standard Armboard
139822 Armboard Pad, Gel, 600[L] x 400[W] x 13mm[H], to suit Wide Armboard



111487

GEL ELBOW PAD

- Positions and cradles the elbow to relieve pressure.
- 85 [dia.] x 38mm[H].

CODE:

111494 Gel Elbow Pad



GEL ULNAR/BRACHIAL PROTECTOR

- Protects the forearm, wraps around the patients arm and secures with hook and loop.
- Helps to prevent ulnar nerve damage when patients arms are tucked to their sides.
- 475[L] x 330[W] x 6mm[H]

CODE:

111490 Gel Ulnar/Brachial Protector



GEL OVAL ULNAR/BRACHIAL PROTECTOR

- Protects elbow and is secured with hook and loop.
- Allows full view of arm.
- 400[L] x 152[W] x 13mm[H].

CODE:

111491 Gel Oval Ulnar/Brachial Protector



GEL ARM AND LEG TOBOGGAN PAD

- Pad to suit surgical arm table guard 111391.

CODE:

111492 Arm and Leg Toboggan Pad, Gel,
430[L] x 190[W] x 10mm[H]



GEL HAND PROTECTOR WITH HOOK AND LOOP

- Encloses hand and secures with hook and loop straps.

CODE:

114660 Hand Protector, Gel, with Hook and Loop,
152[L] x 203mm[W]



Gel Pressure Care Products

GEL LATERAL HEAD PAD WITH CENTRE DISH

- Reduces pressure and protects patient's ear, jaw and side of the face.

CODE:

111615 Gel Lateral Head Pad, Centre Dish, Adult, 254[L] x 228[W] x 25mm[H]



GEL HORSESHOE HEAD PAD

CODE:

111606 Gel Horseshoe Pad, 200[OD] x 45mm[H], Adult
111607 Gel Horseshoe Pad, 83[OD] x 19mm[H], Neo-Natal
111608 Gel Horseshoe Pad, 140[OD] x 32mm[H], Paediatric
111611 Gel Horseshoe Pad, 205[OD] x 76mm[H], Bariatric

*OD = Outer Diameter of Pad.



GEL SUPINE HEAD PAD WITH DISH

- Self centering dish conforms to the patients head while offering support.

CODE:

111614 Gel Supine Head Pad, Adult, 255[L] x 230[W] x 55mm[H]
111672 Gel Supine Head Pad, Paediatric, 177[L] x 152[W] x 50mm[H]
111674 Gel Supine Head Pad, Neo-Natal, 127[L] x 100[W] x 50mm[H]



GEL RICHARD TYPE HEAD REST (SIDE OR SUPINE)

- Improves comfort and stability when patient is lateral or supine.
- 200[L] x 230[W] x 75mm[H]

CODE:

111613 Gel Richard Type Head Rest



GEL HEAD REST PRONE/SUPINE

- Aids in positioning patients with respiration tubes in place.
- Unique design provides superior head control stability.

CODE:

111609 Gel Head Rest Prone/Supine




GEL OPHTHALMIC HEAD PAD

- Contoured shape to support and stabilise the patient's head when in the supine position.
- Designed for general surgery, ophthalmic and oral procedures.
- 254[L] x 305[W] x 98mm[H]

CODE:

116188 Gel Ophthalmic Head Pad



GEL PRONE HEADREST WITH INTERCHANGEABLE BASE

- Including three different foam base heights to provide more optimised positioning.
- Base height: 37mm, 43mm and 62mm.
- 270[L] x 229mm[H]

CODE:

116017 Prone Headrest, Interchangeable Base, Gel




GEL CONTOURED HEAD PAD

- 180[L] x 180[W] x 80mm[H], 32mm[H] contour.

CODE:

111676 Contoured Head Pad, Gel



GEL ROLL COVERS

- Held in place by hook and loop.

CODE:

111715 Gel Roll Cover, Hook and Loop, 430[L] x 300[W] x 13mm[H]
111716 Gel Roll Cover, Hook and Loop, 457[L] x 457[W] x 13mm[H]



GEL DONUT HEAD PAD

CODE:

111603 Gel Donut Head Pad, Adult, 200[OD] x 50mm[H]
111605 Gel Donut Head Pad, Neo-Natal, 90[OD] x 20mm[H]
111604 Gel Donut Head Pad, Paediatric 140[OD] x 32mm[H]
111612 Gel Donut Head Pad, Bariatric, 200[OD] x 75mm[H]



GEL FLAT BOTTOM CHEST ROLLS

- Used to protect the patients chest while in prone position.
- Also used for a variety of positioning needs including oblique and lateral stabilisation.

CODE:

111702 Gel Roll 300[L] x 100[W] x 100mm[H]
111703 Gel Roll 300[L] x 150[W] x 150mm[H]
111704 Gel Roll 400[L] x 100[W] x 100mm[H]
111705 Gel Roll 400[L] x 150[W] x 150mm[H]
111706 Gel Roll 505[L] x 100[W] x 100mm[H]
111707 Gel Roll 505[L] x 150[W] x 150mm[H]



Gel Pressure Care Products

GEL BARIATRIC KNEE CRUTCH PAD

CODE:

- 111499 Gel Bariatric Knee Crutch Pads, Small, 279[L] x 229[W] x 25mm[H]
- 111500 Gel Bariatric Knee Crutch Pads, Medium, 432[L] x 279[W] x 25mm[H]
- 111501 Gel Bariatric Gel Knee Crutch Pads, Large, 508[L] x 508[W] x 25mm[H]



GEL CLAM SHELL HEEL PROTECTOR

- Wraps around the patient's foot to protect the heel.
- Hook and Loop closure.
- Easy to use and clean.

CODE:

- 111594 Gel Clam Shell Heel Protector




GEL HEEL CUP

- Sold individually.

CODE:

- 111595 Gel Heel Cup, Single



GEL KNEE CRUTCH PAD

- Used to protect leg while on knee crutch, Lloyd Davies stirrups and unaffected leg holder.

CODE:

- 111677 Gel Knee Crutch Pad Small 280[L] x 230[W] x 10mm[H]
- 111678 Gel Knee Crutch Pad Medium 430[L] x 280[W] x 6mm[H]
- 111679 Gel Knee Crutch Pad Large 505[L] x 505[W] x 13mm[H]
- 111680 Gel Knee Crutch Pad X-Large 688[L] x 279[W] x 13mm[H]



GEL LITHOTOMY BOOT PAD - BLACK BOOT ALLEN STYLE

- Sold as a pair.
- Gel pads designed to specifically fit Allen stirrup boots.

CODE:

- 111596 Lithotomy Boot Pad, Gel, Small
- 111597 Lithotomy Boot Pad, Gel, Large



GEL LITHOTOMY LEG HOLDER STIRRUP PAD

- For use with lithotomy leg holders, see page 27.

CODE:

- 111600 Leg Holder Stirrup Pad, Gel, 710[L] x 50mm[W], Pair
- 111599 Leg Holder Stirrup Pad, Gel, 710[L] x 50mm[W], Single



GEL LITHOTOMY POLE STRAPS

- For use with lithotomy leg holders, see page 27.

CODE:

- 111601 Lithotomy Pole Straps, Gel, 318[L] x 47[W] x 10mm[H], Pair
- 111602 Lithotomy Pole Straps, Gel, 440[L] x 47[W] x 10mm[H], Pair



GEL FRACTURE TABLE FOOT PAD

CODE:


- 111598 Fracture Table Foot Pad, Gel, 305[L] x 100[W] x 6mm[H]



GEL TRAPEZOID POSITIONER

CODE:

- 111717 Gel Trapezoid Positioner Small, 265[L] x 170[W] x 85mm[H]
- 138396 Gel Trapezoid Positioner, Large, 355[L] x 170[W] x 120mm[H]



CAN'T FIND THE SIZE YOU REQUIRE?

Scan the QR code below or contact our Customer Service Team on 1800 75 93 93 or info@hipac.com.au.



Gel Pressure Care Products

GEL LATERAL POSITIONER

- Consists of a vinyl covered foam support with removable gel sheet.
- Protects arm and shoulder in lateral positioning.

CODE:

111714 Gel Lateral Positioner Pad, 711[L] x 495[W] x 152mm[H]



GEL BARIATRIC LAP BELT PAD

- Used in conjunction with restraint straps.

CODE:

111497 Gel Bariatric Lap Belt Pad, Medium, 915[L] x 150[W] x 6mm[H]
111498 Gel Bariatric Lap Belt Pad, Large, 1220[L] x 175[W] x 6mm[H]



GEL ANDREWS FRAME KNEE PAD

- Use with the Andrews frame to help reduce pressure on the knees.

CODE:

111693 Andrews Frame Knee Pad, Gel, 559[L] x 406[W] x 25mm[H]



GEL PATIENT POSITIONER

CODE:

111708 Gel Positioner 250[L] x 115[W] x 50mm[H]
111709 Gel Positioner 355[L] x 125[W] x 50mm[H]



GEL SANDBAGS

CODE:

111719 Gel Sandbag 150[L] x 115[W] x 50mm[H]
111723 Gel Sandbag 250[L] x 115[W] x 50mm[H]
111724 Gel Sandbag 300[L] x 115[W] x 50mm[H]
111725 Gel Sandbag 355[L] x 125[W] x 50mm[H]



GEL BASSINET PAD

CODE:

111689 Bassinet Pad, Gel,
660[L] x 305[W] x
13mm[H]



GEL ISOLETTE PAD

CODE:

111690 Isolette Pad, Gel,
660[L] x 405[W] x
13mm[H]



GEL PAEDIATRIC PAD

CODE:

111688 Paediatric Pad, Gel, 280[L] x
230[W] x 10mm[H], Small
111687 Paediatric Pad, Gel, 432[L] x
279[W] x 6mm[H], Medium



GEL MAYO STAND COVER

CODE:

111695 Mayo Stand Cover,
Gel, 530[L] x 470[W] x
16mm[H]



GEL EPIDURAL POSITIONING DEVICE ARM REST PAD

CODE:

111694 Epidural Device Arm
Rest Pad, Gel, 330[L] x
127[W] x 6mm[H]



GEL CHEST POSITIONERS

- Contoured shape conforms to the patient's chest and shoulder.
- Sold as a single unit or as a set.
- 711[L] x 203[W] x 76mm[H] each.

CODE:

111710 Chest Positioners, Gel, Left
111711 Chest Positioners, Gel, Right
111712 Chest Positioners, Gel, Set
[Left and Right]



Gel Pressure Care Products

GEL SUMMERSVILLE PAD

CODE:

111700 Summersville Pad, Gel, 304[L] x 304[W] x 6mm[H]



GEL TABLE PEG BOARD PAD SET

- Gel pad set to suit peg board system.
- One pad: 1200[L] x 510[W] x 12.5mm[H].
- Two pads: 200[L] x 200[W] x 10mm[H].
- Two pads: 190mm[L] x 130mm[W] x 10mm[H].

CODE:

114114 Table Peg Board Pad Set, Gel



GEL PAD TO SUIT WILSON FRAME

- Sold individually.
- 737[L] x 254[W] x 13mm[H].

CODE:

111691 Gel Pad to suit Wilson Frame



GEL TABLE HIP PAD

CODE:

111701 Table Hip Pad, Gel, 505[L] x 508[W] x 13mm[H]



GEL SURGICAL ARM TABLE PAD

CODE:

119603 Surgical Arm Table Pad, Gel, 760[L] x 440[W] x 13mm[H]
139821 Surgical Arm Table Pad, Gel, 800[L] x 400[W] x 13mm[H]



GEL HOURGLASS ARM TABLE PAD

CODE:

139823 Hourglass Arm Table Pad, Gel, 800[L] x 540[W] x 13mm[H]



GEL SACRAL PAD

- Gel pad to suit sacral rest.
- With perineal post hole.

CODE:

116968 Sacral Pad, Gel, 228[W] x 178[L] x 13mm[H]



GEL RESTRAINT STRAP WITH CLIPS AND BUCKLE

- 2400mm long restraint strap including clips and buckle, with gel lap pad insert.
- Pad dimensions: 914[L] x 152[W] x 6mm[H].

CODE:

137909 Restraint Strap with Clips and Buckle, Includes Gel Insert



GEL DOMED POSITIONER

CODE:

138704 Domed Positioner, Gel, 370[L] x 160[W] x 70mm[H]



GEL SAMPLE PAD

CODE:

111697 Sample Pad, Gel, 125[L] x 100mm[W]



GEL REPAIR PATCH KIT

- Easily repairs cuts, punctures and cracks.
- Four patches per kit.

CODE:

111696 Gel Repair Patch Kit, 100 x 100mm

