



Mayo instrument stand, designed for use in the Operating Room, offers all the quality and safety standards necessary in this environment. Made entirely out of 18/10 AISI 304 Stainless Steel [X5CrNi1810], it ensures,

- High level of resistance to corrosion
- High level of bacterial removal
- Low level of bacterial retention
- Excellent durability

Features

- Solid 3 spoke rectangular section tubular base, with safety tips equipped with 3 twin castors, 75mm, with threadguard, in non-marking polymer, of which 2 are breaking and 1 is anti-static.
- Removable deep-drawn table top with closed, rounded corners. Degree of freedom $+180^\circ / -180^\circ$, degree of operating safety $+90^\circ / -90^\circ$. Rotation system equipped with safety lock.
- Lift system with double acting high precision hydraulic pump, with sealed circuit, pedal-operate.
- Lift range from 910mm to 1270mm.
- Pedal in ergonomic position placed inside the structure.
- The materials, components and assembly systems allow disinfection and cleaning in trolley washer disinfectant.
- Maximum load capacity 30kg, evenly distributed on the table top.

CODE

142252

Table top dimension 750mm x 550mm x 30mm
Outer dimensions 780mm x 580mm x 910mm

142487

Table top dimensions 650mm x 530mm x 20mm
Outer dimensions 780mm x 580mm x 910mm

ABOUT HIPAC

We are an Australian healthcare company that guarantees to stand behind all our medical products, both manufactured and imported, because we're here for the life of our products and clients.

PO Box 797, 36 Long St. Goulburn NSW 2580, Australia
P. +61 2 4823 0000 E. info@hipac.com.au

SYDNEY | MELBOURNE | BRISBANE | ADELAIDE | PERTH

hipac.com.au

ABN 27 600 353 688