

Hydraulic Surgeon Stool**Datasheet****Introduction**

The hydraulic surgeon stool is suitable for a broad range of surgical specialties. The advanced designs provide many features and advantages. The hydraulic foot pedals allows for easy and hands-free height adjustment. The backrest and armrests are also adjustable.

Features

- Adjustable contoured backrest
- Dual-density foam seat cushion
- Foot-controlled, hydraulic, three spindled base
- Front-end locking system
- ABS thermoformed plastic hood
- 360° rotate armrest pad, 180° armrest support
- Height adjustment: 530mm-710mm
- Weight Capacity: 136kg



Water Drop



Wedge



Saddle

Foot-controlled hydraulic [ascend]



Front-end locking system

Foot-controlled hydraulic [decline]