

**LUCINA 4**

**PRODUCT BROCHURE**



## Lucina 4 Overview

Simple LED strobe-guided positioning allows caregivers to focus on the patient.

### Remote Controlled Positioning

Featuring Skytron exclusive strobe-guided technology, the Lucina 4 positioning wand ensures that the light and caregiver remain focused on the patient. Guide the light beam to the targeted location with a simple click of a button. Available in configurations ranging from one to four lights, the Lucina 4 wand controls positioning and on/off of all lights simultaneously.



### Pleasing Room Aesthetics

Lucina 4 lights are mounted flush with the ceiling, ensuring that each clinical space remains comfortable and less intimidating for the patient. In areas with low ceiling heights or other mounting restrictions, the Lucina 4 integrates seamlessly and is less invasive than traditional exam lighting arm configurations. Works well as a diagnostic/exam light in LDR, ICU, Hybrid Operating Rooms, Cath Labs, Exam Rooms and more.

### Energy Efficient, Long Lasting LEDs

Combined with Skytron's segmented reflector design the Lucina 4 produces brilliant, colour correct lighting for a variety of clinical settings.

## Lucina 4 Features

### Features

- Unique LED strobe-guided positioning wand
- Angled up to 47° for superior illumination
- Long lasting LEDs
- Fully recessed 24" square fixture
- Wall control with ON/OFF and intensity touchpad and (back-up) positioning buttons
- Available in one, two, three and four fixture configurations



### On-Demand Accuracy in any Setting

The Lucina 4 is perfect for bringing simple operation and flexible positioning to any clinical environment, helping keep healthcare spaces open and free of ceiling mounted equipment conflicts.

#### Ideal for:

- Diagnostic and Exam Light
- Labour, Delivery, Recovery and Postpartum Suites
- Intensive Care Units
- Hybrid Operating Rooms
- Cath Labs

Lucina 4 is a diagnostic/examination light and is not intended to be used as a primary surgical luminaire. Hipac recommends a surgical light (fixed or mobile) be included in clinical settings where surgical illumination may be required.

## **ABOUT HIPAC**

We are an Australian healthcare company that guarantees to stand behind all our medical products, both manufactured and imported, because we're here for the life of our products and clients.

PO Box 797, 36 Long St. Goulburn NSW 2580, Australia  
P. +61 2 4823 0000 or 1800 75 93 93  
E. [info@hipac.com.au](mailto:info@hipac.com.au)

SYDNEY | MELBOURNE | BRISBANE | ADELAIDE | PERTH

**[hipac.com.au](http://hipac.com.au)**

ABN 27 600 353 688