

SPECTRA EXAMINATION LIGHTS

PRODUCT BROCHURE



Reliable, Versatile and Cost Saving

Featuring seven LEDs with unique vertically segmented reflectors, the Spectra exam light delivers crisp illumination with reduced shadows. In addition to the brilliant spot, caregivers will enjoy the flexible positioning and drift-free stability- bringing confidence and fewer hassles to any patient exam.

The durable construction and sealed, easy to clean lighthead make the Spectra exam light ready for any healthcare environment. Bring quality illumination to the Emergency, Critical Care, Labour and Delivery, imaging department and more



Features

- Vertically segmented reflectors deliver a clear spot with fewer shadows.
- Brilliant, white light is 51, 000 Lux at one metre
- 355° lighthead pitch extends upward and downward positioning capabilities
- Cool, colour-correct light provides CRI of 96
- Long lasting LEDs
- Low-profile design for simple transport and storage
- · Optional sterilisable handle available.



Spectra Ceiling Mounted Examination Light

- Ceiling or boom mounted design
- 71" reach
- 360° arm rotation
- 355° lighthead rotation





Spectra Wall Mounted Examination Light

- Wall mounted design
- 75.5" reach
- 180° arm rotation
- 355° lighthead rotation



Spectra Mobile Stand Examination Light

- Mobile Stand
- 38" reach
- 50° upward lighthead travel
- 30° downward lighthead travel



ABOUT HIPAC

We are an Australian healthcare company that guarantees to stand behind all our medical products, both manufactured and imported, because we're here for the life of our products and clients.

PO Box 797, 36 Long St. Goulburn NSW 2580, Australia P. +61 2 4823 0000 or 1800 75 93 93 E. info@hipac.com.au

SYDNEY | MELBOURNE | BRISBANE | ADELAIDE | PERTH

hipac.com.au

ABN 27 600 353 688