

STORAGE AND WORKSTATION SOLUTIONS

PRODUCT CATALOGUE



Select Series Anaesthetic Cart



Select Series Anaesthetic Cart and Accessory Package

CODE:	COLOUR:	LOCKING MECHANISM:
118837	Slate Blue	Key Lock

- 1 x 117788 Select Series Anaesthetic Cart
- 1 x 113137 Select Series Trellis
- 1 x 113132 Select Series Tilt Bins Five Units
- 1 x 113131 Select Series Tilt Bins Four Units
- 1 x 113286 Tubular Wrap Around Railing
- 1 x 113212 Glove Box Holder
- 1 x 113298 Tape Dispenser with Brackets

SIZE:

Cart: 812 (L) x 635(W) x 1193(H) Drawer: 3 x 76mm, 2 x 152mm, 1 x 228mm

- 1 x 113236 Small Draw Divider Pack Includes Two x 76mm(D) Trays, Four Long Dividers, Eight Short Dividers
- 1 x 113240 Large Draw Divider Pack Includes Two x 127mm(D) Trays, Two Long Dividers, Four Short Dividers
- 1 x 113327 Waste Container w/Lift Top
- 1 x 113321 Wire Basket 305(L) x 305(W) x 203mm(D)
- 1 x 113223 76mm(D) Tray with 8 Flex
 Dividers and 14 Ampule Holders

Select Series Anaesthetic Cart (No Accessories)

CODE:	COLOUR:	LOCKING MECHANISM:	SIZ
117788	Slate Blue	Key Lock	Ca

SIZE:

Cart: 812(L) x 635(W) x 1193mm(H) Drawer: 3 x 76mm, 2 x 152mm, 1 x 228mm

Select Series Anaesthetic Cart Accessories

TUBULAR WRAP AROUND RAILING

CODE: 113286



LARGE DRAWER DIVIDER PACK CODE: 113240 SIZE: 419(L) x 279(W) x 127mm(D) Includes Two Trays, Two Long Dividers, Four Short Dividers



SMALL DRAWER DIVIDER PACK

CODE: 113236 SIZE: 419(L) x 279(W)x 76mm(D) Includes Two Trays, Four Long Dividers, Eight Short Dividers



SELECT SERIES TILT BINS

CODE: 113131 - Four Units 113132 - Five Units



TRAY WITH FLEX DIVIDERS

CODE: 113223 SIZE: 76m(D) Tray with 8 Flex Dividers and 14 Ampule Holders



WIRE BASKET CODE: 113321 SIZE: 305(L) x 305(W)x 203mm(D)



We're here for life

WASTE CONTAINER WITH LIFT TOP





SELECT SERIES TRELLIS CODE: 113137



TAPE DISPENSER WITH BRACKETS CODE: 113298



GLOVE BOX HOLDER CODE: 113212



Select Series Emergency Cart



Select Series Anaesthetic Cart and Accessory Package

CODE: COLOUR: 118836 Red

LOCKING MECHANISM:

Breakaway Lock

- 1 x 117787 Select Series Emergency Cart
- 1 x 113214 Height Adjustable IV Pole
- 1 x 113250 Oxygen Tank Brackets
- 1 x 112464 Select Series Cardiac Board

SIZE:

Cart: 812(L) x 635(W) x 1193mm(H) Drawer: 3 x 76mm, 2 x 152mm, 1 x 228mm

- 1 x 113236 Small Draw Divider Pack Includes Two x 76mm(D) Trays, Four Long Dividers, Eight Short Dividers
- 1 x 113240 Large Draw Divider Pack Includes One x 127mm(D) Trays, Two Long Dividers, Four Short Dividers
- 1 x 113261 Plastic Breakaway Plastic Seals 100/pk

Select Series Anaesthetic Cart (No Accessories)

CODE:	COLOUR:	LOCKING MECHANISM:
117787	Red	Breakaway Lock

SIZE:

Cart: 812(L) x 635(W) x 1193mm(H) Drawer: 3 x 76mm, 2 x 152mm, 1 x 228mm

Select Series Emergency Cart Accessories

TUBULAR WRAP AROUND RAILING CODE: 113286



LARGE DRAWER DIVIDER PACK CODE: 113240

SIZE: 419(L) x 279(W) x 127mm(D) Includes Two Trays, Two Long Dividers, Four Short Dividers



SMALL DRAWER DIVIDER PACK

CODE: 113236 SIZE: 419(L) x 279(W) x 76mm(D) Includes Two Trays, Four Long Dividers, Eight Short Dividers

BREAKAWAY PLASTIC SEALS

CODE: 113261 100 pack



HEIGHT ADJUSTABLE IV POLE CODE: 113214

WIRE BASKET CODE: 1113321

SI7F:



• We're here for life

OXYGEN CYLINDER BRACKET CODE: 113250



SELECT SERIES SHARPS BRACKET

CODE: 113128 Bracket only- does not include sharps bin



SELECT SERIES CARDIAC BOARD CODE: 112468



Select Series Paediatric Cart



Select Series Paediatric Cart

CODE:	COLOUR:	LOCKING MECHANISM:
117789	Red	Beige, Pink/Red, Lavender,
		Yellow, White, Crash Cart
		Blue, Orange, Lime Green
		and Beige

NOTE: See Pages Four and Six for Accessory Options

- Twin poly handles
- Durable sealed work surface
- Full extension ball bearing drawer
- Outrigger stabilising frame with bump
- Four premium 125mm Castors; two locking, one swivel, one directional
- Full length breakaway lock

Cart: 812[L] x 635[W] x 1193mm[H] Drawer: 8 x 76mm, 1 x 152mm

- Single pullout side shelf
- Colour coded drawers

SIZE:



-We're here for life



Features:

- Key lock on all drawers
- Four swivel castors, two lockable
- Durable sealed work surface

Select Series Mini Value Cart

CODE:	COLOUR:	LOCKING MECHANISM:
111990	Taupe	Key Lock

SIZE:

Cart: 457(L) x 457(W) x 863mm(H) Drawer: 3 x 76mm, 1 x 152mm, 1 x 229mm

ABOUT HIPAC

We are an Australian healthcare company that guarantees to stand behind all of our medical products, both manufactured and imported, because we're here for the life of our products and clients.

PO Box 797, 36 Long Street. Goulburn NSW 2580, Australia P. +61 4823 0000 or 1800 75 93 93 E. info@hipac.com.au

SYDNEY | MELBOURNE | BRISBANE | ADELAIDE | PERTH

hipac.com.au

ABN 27 600 353 688