

3503 EZ Slide Operating Table

The Standard in Full Body Imaging and Patient Positioning



The 3503 EZ Slide features a convenient top slide and the highest table height in our portfolio. With 530mm of top slide, imaging and surgical access are easy, allowing patient movement through the imaging window with the click of a button. Its 1092mm table height provides comfortable working access for a broad range of surgical specialties.

Hipac's extensive range of table accessories and traction frames are compatible with our full range of operating tables.

- comfortable surgical access
- powerful articulation capacity
- one-touch imaging access
- the 530mm top slide allows greater patient movement through the imaging window.
- 320kg weight capacity with full top slide & 272kg articulated weight capacity with full top slide
- removable head and leg sections for quick accessory attachment and positioning
- 20° lateral tilt, 30° Trendelenburg and 30° reverse Trendelenburg
- slim profile base that stays clear of mobile imaging equipment, providing hassle free C-arm access
- quick integration with accessories such as our Orthopaedic Traction Frame
- supplied with side rails of choice, so that existing accessories will be immediately usable, avoiding unnecessary and often expensive replacement accessories

One-Touch Pendant Control

- simple one-touch operation
- backlit control for high visibility in low lighting
- audible alarm for kidney bridge, back up articulation and top slide obstruction
- with PSR (Patient Safe Return)
- recessed buttons prevent accidental activation
- durable and water resistant
- can be used in conjunction with optional three function foot control




Back-up controls for emergency situations



Battery life/diagnostic indicator



Easy to use clip style power cord



Optional 40" carbon fiber leg section for lower body imaging



Back lit one-touch positioning functions include Trendelenburg, tilt, and more

Positioning Guide



Upper Body Imaging



Lower Body Imaging



Lap Nissen



Spinal



Ophthalmic/ENT



Nephrectomy



Neuro (Neck)



Neuro (lumbar)



Shoulder Arthroscopy



Urology



Lower Extremity



Bariatric


3503 EZ SLIDE OPERATING TABLE



635mm Minimum
Table Height



1092mm Maximum
Table Height



530mm Total
Slide



40° Flex



90° Reflex



98° Leg Down



90° Back Up



40° Back Down



30° Trendelenburg



30° Reverse
Trendelenburg



20° Tilt - Left



20° Tilt - Right

ABOUT HIPAC

Operating at the forefront of healthcare, we strive to redefine what's possible. We provide superior products, design expertise and total solutions to help create environments where nothing stands in the way of patient care.

Contact information



1800 75 93 93



info@hipac.com.au



www.hipac.com.au



We're here for life