



SPECTRA[®]

Our premium Truly Cordless[™]
surgical headlight system.





Our Features



Modular Bezels

Upgrade or change out your bezel with ease with our bezel port.




Convertible Top

Each Spectra comes with both a hard top and a soft top to choose from for optimal comfort.



Comfort Headband

Comfort is a necessity. With a patented padding system, get the comfort needed and forget it's even there.



Brightness

Choose a light that can power a beam bright enough to see all you need when it matters.



Adjustable

Focus beam in your desired work area.



Truly Cordless

Cut the cord and heavy belt battery packs for good. With our lightweight Truly Cordless™ lithium batteries, cords will be a thing of the past.

Our Headlight Kits

Our premium Truly Cordless™ surgical headlight system.



MG8 Bezel



MG3 Bezel



MG6 Bezel

	Product SKU	Kit Type	Headlight	45w Wall Plug / USB	Padding Conversion Kit	Carrying Case	Batteries	4-Bay Charger
MG3	MG3-SK	Standard Kit	✓	✓	✓	✓	x2	
E-Class Clinical	MG3-HK	Hospital Kit	✓	✓	✓	✓	x6	✓
MG6	MG6-SK	Standard Kit	✓	✓	✓	✓	x2	
Crossover Clinical & O.R.	MG6-HK	Hospital Kit	✓	✓	✓	✓	x6	✓
MG8	MG8-SK	Standard Kit	✓	✓	✓	✓	x2	
Premium Operating Room	MG8-HK	Hospital Kit	✓	✓	✓	✓	x6	✓

Our Bezels



MG3

E-Class
Clinical




Battery Life
10 hours



Brightness
100,000 LUX



Spot Size
10-100 mm



MG6

Crossover
Clinical & O.R.




Battery Life
8 hours



Brightness
200,000 LUX



Spot Size
70-90 mm



MG8


Premium
Operating Room




Battery Life
5 hours



Brightness
300,000 LUX



Spot Size
100-200 mm



Accessories

Spare Batteries & Battery Chargers



Description

Battery Set (2 pack)
 Battery Charger 4-Bay Desktop Console Set (US)
 Battery Charger 4-Bay Desktop Console Set (International)

Product SKU

M-90007
 M-90206
 M-90207

Carrying Case & Bag



Description

Hard Carrying Case
 Carrying Bag - Rugged
 Carrying Bag - Slim Drawstring

Product SKU

M-90912
 M-90910
 M-90907

Plugs & Adapters



Description

Battery Charger Wall Set (US)
 Battery Charger Wall Set (International)

Product SKU

M-90102
 M-90106

Magnifiers lenses - MG3 (small)



Diopter	Magnification	Increase in eyesight	Optical power	Recommended working distance (focal length)	Product SKU
+2	+0.40	40%	1.4x	250 mm	M-90301
+2.5	+0.50	50%	1.5x	235 mm	M-90302
+3	+0.60	60%	1.6x	220 mm	M-90303
+3.5	+0.70	70%	1.7x	210 mm	M-90304
+4	+0.80	80%	1.8x	200 mm	M-90305
+4.5	+0.90	90%	1.9x	195 mm	M-90306

Magnifiers lenses MG6 / MG8 (large)



Diopter	Magnification	Increase in eyesight	Optical power	Recommended working distance (focal length)	Product SKU
+2	+0.40	40%	1.4x	250 mm	M-90321
+2.5	+0.50	50%	1.5x	235 mm	M-90322
+3	+0.60	60%	1.6x	220 mm	M-90323
+3.5	+0.70	70%	1.7x	210 mm	M-90324
+4	+0.80	80%	1.8x	200 mm	M-90325
+4.5	+0.90	90%	1.9x	195 mm	M-90326

Accessories



Bezel Systems

Description	Product SKU
MG3 Bezel System	M-92431
MG6 Bezel System	M-92461
MG8 Bezel System	M-92471



Padding Sets

Description	Product SKU
Padding Set (4 pc.) - Snap-Around	M-91401
Headstrap Top Kit (Hard-Top & Soft-Top)	M-91405
Padding Rear (1 pc.) - Snap AirPad Assembly	M-91411
Headstrap	M-92401



Face Shield

Description	Product SKU
Face Shield & Mounting Kit	M-91420
Face Shield & Mounting Kit (with extra Shield)	M-91420-2
Face Shield - Standalone Headband for Loupes & Loupe Lights	M-91422
Face Shield - Standalone Headband for Loupes & Loupe Lights (with extra Shield)	M-91422-2
Padding Set Replacement for Stand-alone Face Shield	M-91425
Mounting Knob Set Replacement for MedLED Headlight	M-91426
Shield (Rigid) Replacement for MedLED Headlight	M-91427
Shield (Rigid) Replacement for Loupes & Loupe Lights	M-91428
Shield (flex) Replacement for Loupes & Loupe Lights	M-91428



Amber Filters

Description	Product SKU
Amber Filter - (for MG3)	M-91101
Amber Filter - (for MG6)	M-91103



Storage Options

Description	Product SKU
4-bay Wall Mount Storage Unit	M-90711

SPECTRA[®]



CORPORATE HEADQUARTERS
825 NE 25th Avenue
Portland, OR 97232 USA

Office: +01 (503) 517-8001
Toll Free: +01 (866) 445-4923
Fax: +01 (503) 517-8002
E-Mail: info@rpmed.com

www.rpmed.com



www.medled.com

ISO 13485:2016 | FDA Registered