Hipac
We're here for life

# Life-Line Key Wear

## Inspection Guide

## KEY WEAR REFERENCE

Compare the condition of Life-Line keys against the examples below.



**NEW KEY**



**PART WORN KEY**

(Order replacement key)




**WORN KEY**

(Replace key)

### IMPORTANT NOTE:


The Life-Line Keys are designed to wear out before the lock.
Key wear should be constantly monitored before failure occurs.

## REORDER REPLACEMENT KEYS




**120337**

LifeLine Butterfly Key




**137322**

Small Version LifeLine Butterfly Key


### ABOUT HIPAC

We are an Australian healthcare company that guarantees to stand behind all our medical products, both manufactured and imported, because we're here for the life of our products and clients.

PO Box 797, 36 Long St. Goulburn NSW 2580, Australia
P. +61 2 4823 0000 E. info@hipac.com.au

SYDNEY | MELBOURNE | BRISBANE | ADELAIDE | PERTH

**hipac.com.au**

ABN 27 600 353 688



Hipac
We're here for life