

SKU: 148171

Hipac SureLift Stirrups

Hipac's SureLift Stirrup system safely positions and supports patients' legs during lithotomy procedures, improving workflow efficiency, clinician ergonomics, and patient safety.

Designed for use in operating theatres, SureLift Stirrups enable smooth intraoperative adjustment of lithotomy and abduction positions while maintaining sterility, stability, and procedural confidence.

Suitable for a wide range of procedures including:

- Gynae/Uro
- Robotic Surgery
- General Surgery
- Cystoscopy
- D&C Procedures
- Laparoscopic Surgery
- Colorectal Procedures
- Oncology Procedures
- LEEP Procedures
- TAVH / LAVH Procedures

Hipac
We're here for life

Key Features

The anatomically contoured boot design distributes pressure evenly across the lower leg, helping reduce localised pressure points while securely supporting the limb. The lateral support fin assists in maintaining proper knee alignment during higher degrees of lithotomy, helping minimise strain through the groin and lower back. The reusable padded boot fully supports the foot, ankle and calf to enhance comfort and positioning stability during procedures.

Assisted Lift Technology:
Helps reduce the effort required to position and reposition the patient's leg, supporting safer manual handling practices.



Floating Boot Design:
Self-adjusting support helps minimise calf pressure across varying lithotomy positions.

Single-Handed Adjustment Handle:
Allows quick and controlled intraoperative positioning adjustments while maintaining the sterile field.



Multi-Axis Positioning:
Enables simultaneous lithotomy and abduction adjustment for efficient and precise patient positioning.

Position Indicators:
Visual indicators assist with repeatable and accurate positioning.



Enhanced Patient Support:
Boot-style support helps distribute pressure more evenly than traditional stirrup designs.

Hipac

Contact Information

1800 75 93 93 | info@hipac.com.au | www.hipac.com.au


Specifications

Technical Details	
Lithotomy Range	+84° to -33°
Abduction Range	+25° to -9°
Patient Weight Limit	227kg



Accessories


Stirrup Clamps	
A Rail	SKU: 113091
B Rail	SKU: 113102
D Rail	SKU: 113111
E Rail	SKU: 113119



Lithotomy Storage Cart
SKU: 119678



Universal Bar Clamps	
A Rail	SKU: 115167
B Rail	SKU: 115699
D Rail	SKU: 115700
E Rail	SKU: 115701



Universal Stirrup Trolley
SKU: 120728



To view our complete list of products, scan the QR code or visit www.hipac.com.au

Hipac

Contact Information

1800 75 93 93 | info@hipac.com.au | www.hipac.com.au